

# Health

## MEDIFLEX



### DESCRIPTION

CORE	High - Resilience polyurethane foam 38 or 43 kg/cubic meter, cut in a waffle-pattern on the «sleeper» side to prevent bed-sores, multi articulated on the other side.
FABRIC	Polyester knitted cover polyurethane coating, permeable, green color (150 g/sqm)
FINISH	Removable cover, zippered on the bottom edge on 1 length and 1 width, with an identification label on the «sleeper» side.
WASHING	Machine-washable at 90°C, tumble-drying is acceptable up to 90°C, use soapy water, and/or very diluted bleach-water. Do NOT use ether or other softeners
TREATMENT	Fabric treated with anti-bacterial, anti-fungal, solution, permeable (400g/sqm/24H)
CONFORMITY	Mattress conforms NF EN 597, parts 1 and 2
COMFORT INDEX	Compression resistance at 40% (kPa): 3.9 (38 kg/cubic meter) or 3.9 (43 kg/cubic meter)
CLIENTELE	Mattress for patients with a Norton score $\geq 14$ (low-risk), or stage 1 bed-sores
GUARANTEE	Against premature sagging
OPTION	Bevelled corners

### LABELS

